

Table of Contents

Chapter No.	Chapter Name	Page No.
C 1	Hand & Foot Plates	5
C 2	Locking Distal Radius Plates	7-9
C 3	2.7 mm Locking Plates	11-12
C 4	3.5 mm Locking Plates	15-24
C 5	4.5/5 mm Locking Plates	27-32
C 6	DHS-DCS & Angle Plates	35-37
C 7	Bone Screws	39-50
C 8	CranioMaxilloFacial	53-60
C 9	Interlocking Nails	62-73
C 10	Hip Prosthesis	76-77
C 11	Spine System	80-82
C 12	Pins-Wires-Staples	84-92
C 13	External Fixators	94-105

C1 Locking Hand & Foot Plates

1. 2/2.4/2.7 mm Locking Compression Plate
2. 2/2.4/2.7 mm Locking Adaption Plate
3. 2/2.4/2.7 mm Locking T Adaption Plate 3 Head Holes
4. 2/2.4/2.7 mm Locking T Adaption Plate 2 Head Holes
5. 2/2.4/2.7 mm Locking Y Adaption Plate
6. 2/2.4/2.7 mm Locking Condylar Plate
7. 2/2.4/2.7 mm Locking L Plate
8. 2/2.4/2.7 mm Locking L Plate Angle

1.1 2/2.4/2.7 mm Locking Compression Plate



Dia (mm)	2 2.4 2.7
Holes	4, 5, 6, 7, 8, 10
Material	Titanium

1.2 2/2.4/2.7 mm Locking Adaption Plate



Dia (mm)	2 2.4 2.7
Holes	4, 5, 6, 7, 8, 10
Material	Titanium

1.3 2/2.4/2.7 mm Locking T Adaption Plate 3 Head Holes



Dia (mm)	2 2.4 2.7
Holes	4, 5, 6, 7, 8, 10
Material	Titanium

1.4 2/2.4/2.7 mm Locking T Adaption Plate 2 Head Holes



Dia (mm)	2 2.4 2.7
Holes	4, 5, 6, 7, 8, 10
Material	Titanium

1.5 2/2.4/2.7 mm Locking Y Adaption Plate



Dia (mm)	2 2.4 2.7
Holes	4, 5, 6, 7, 8, 10
Material	Titanium

1.6 2/2.4/2.7 mm Locking Condylar Plate



Dia (mm)	2 2.4 2.7
Holes	4, 5, 6, 7, 8, 10
Material	Titanium

1.7 2/2.4/2.7 mm Locking L Plate



Dia (mm)	2 2.4 2.7
Holes	3, 4
Direction	Left Right
Material	Titanium

1.8 2/2.4/2.7 mm Locking L Plate Angle



Dia (mm)	2 2.4 2.7
Holes	3, 4
Direction	Left Right
Material	Titanium

C2 Locking Distal Radius Plates

1. 2.7 mm Variable Angle Locking Two Column Volar Distal Radius Plate
2. 2.7 mm Locking Dorsal Distal Radius Straight Plate
3. 2.7 mm Locking Volar Column Distal Radius Plate 9 Head Holes
4. 2.7 mm Locking Volar Column Distal Radius Plate 8 Head Holes
5. 2.7 mm Locking Dorsal Distal Radius T Plate 3 Head Holes
6. 2.7 mm Locking Volar Long Distal Radius Plate 4 Head Holes
7. 2.7 mm Locking Volar Extra-Articular Distal Radius Plate 4 Head Holes
8. 2.7 mm Locking Volar Extra-Articular Distal Radius Plate 5 Head Holes
9. 2.7 mm Locking Volar Juxta-Articular Distal Radius Plate 5 Head Holes
10. 2.7 mm Locking Volar Distal Radius Buttress Plate 5 Head Holes
11. 2.7 mm Locking Dorsal Distal Radius L Plate 2 Head Holes
12. 2.7 mm Locking Dorsal Distal Radius L Plate 3 Head Holes
13. 2.7 mm Locking Dorsal Distal Radius L Plate Oblique
14. 2.7 mm Locking Proximal Radius Radial Head Rim Plate
15. 2.7 mm Locking Proximal Radius Radial Head Neck Plate

2.1 2.7 mm Variable Angle Locking Two Column Volar Distal Radius Plate



Dia (mm)	2.7
Holes	3, 4, 5
Direction	Left Right
Material	Titanium

2.2 2.7mm Locking Dorsal Distal Radius Straight Plate



Dia (mm)	2.7
Holes	5, 6
Material	Titanium

2.3 2.7 mm Locking Volar Column Distal Radius Plate 9 Head Holes



Dia (mm)	2.7
Holes	3, 4, 5
Direction	Left Right
Material	Titanium

2.4 2.7 mm Locking Volar Column Distal Radius Plate 8 Head Holes



Dia (mm)	2.7
Holes	3, 4, 5
Direction	Left Right
Material	Titanium

2.5 2.7 mm Locking Dorsal Distal Radius T Plate 3 Head Holes



Dia (mm)	2.7
Holes	3, 4, 5, 6
Material	Titanium

2.6 2.7 mm Locking Volar Long Distal Radius Plate 4 Head Holes



Dia (mm)	2.7
Holes	8, 10, 12
Material	Titanium

2.7 2.7 mm Locking Volar Extra-Articular Distal Radius Plate 4 Head Holes



Dia (mm)	2.7
Holes	3, 4, 5
Direction	Left Right
Material	Titanium

2.8 2.7 mm Locking Volar Extra-Articular Distal Radius Plate 5 Head Holes



Dia (mm)	2.7
Holes	3, 4, 5
Direction	Left Right
Material	Titanium

2.9 2.7 mm Locking Volar Juxta-Articular Distal Radius Plate 5 Head Holes



Dia (mm)	2.7
Holes	3, 4, 5
Direction	Left Right
Material	Titanium

2.10 2.7 mm Locking Volar Distal Radius Buttress Plate 5 Head Holes



Dia (mm)	2.7
Holes	3
Direction	Left Right
Material	Titanium

2.11 2.7 mm Locking Dorsal Distal Radius L Plate 2 Head Holes



Dia (mm)	2.7
Holes	3, 4, 5
Direction	Left Right
Material	Titanium

2.12 2.7 mm Locking Dorsal Distal Radius L Plate 3 Head Holes



Dia (mm)	2.7
Holes	3, 4, 5
Direction	Left Right
Material	Titanium

2.13 2.7 mm Locking Dorsal Distal Radius L Plate Oblique



Dia (mm)	2.7
Holes	3, 4
Direction	Left Right
Material	Titanium

2.14 2.7 mm Locking Proximal Radius Radial Head Rim Plate



Dia (mm)	2.7
Holes	2, 3, 4
Direction	Left Right
Material	Titanium

2.15 2.7 mm Locking Proximal Radius Radial Head Neck Plate



Dia (mm)	2.7
Holes	2, 3, 4
Material	Titanium

C3 2.7mm Locking Plates

1. 2.7 mm Locking Compression Plate
2. 2.7 mm Locking Reconstruction Plate Straight
3. 2.7 mm Locking Reconstruction Plate Curved
4. 2.7 mm Locking Reconstruction Plates Straight Low Profile
5. 2.7 mm mm Locking Reconstruction Plates Curved Low Profile
6. 2.7 mm Locking T Plate 2 Head Holes
7. 2.7 mm Locking T Plate 3 Head Holes
8. 2.7 mm Locking T Plate Angle 3 Head Holes
9. 2.7 mm One Third Tubular Plate
10. 2.7 mm Locking Pediatric Hip Plates

3.1 2.7 mm Locking Compression Plate



Dia (mm)	2.7 mm
Holes	4, 5, 6, 7, 8, 10, 12, 14
Material	SS316L, Titanium

3.2 2.7 mm Locking Reconstruction Plate Straight



Dia (mm)	2.7 mm
Holes	4, 6, 8, 10, 12, 14, 16
Material	SS316L, Titanium

3.3 2.7 mm Locking Reconstruction Plate Curved



Dia (mm)	2.7 mm
Holes	4, 6, 8, 10, 12, 14, 16
Material	SS316L, Titanium

3.4 2.7 mm Locking Reconstruction Plates Straight Low Profile



Dia (mm)	2.7 mm
Holes	4, 6, 8, 10, 12, 14, 16
Material	SS316L, Titanium

3.5 2.7 mm Locking Reconstruction Plates Curved Low Profile



Dia (mm)	2.7 mm
Holes	4, 6, 8, 10, 12, 14, 16
Material	SS316L, Titanium

3.6 2.7 mm Locking T Plate 2 Head Holes



Dia (mm)	2.7 mm
Holes	3, 4, 5
Material	SS316L, Titanium

3.7 2.7 mm Locking T Plate 3 Head Holes



Dia (mm)	2.7 mm
Holes	3, 4, 5
Material	SS316L, Titanium

3.8 2.7 mm Locking T Plate Angle 3 Head Holes



Dia (mm)	2.7 mm
Holes	3, 4, 5
Direction	Left Right
Material	SS316L, Titanium

3.9 2.7 mm Locking One Third Tubular Plate



Dia (mm)	2.7 mm
Holes	4, 6, 8, 10, 12, 14, 16, 18
Material	SS316L, Titanium

3.10 2.7 mm Locking Pediatric Hip Plate



Size (mm)	2.7
Holes	100° (3 holes) 110° (3 holes) 120° (4 holes) 130° (3, 5, 7 holes) 150° (3 holes)
Angle	100°, 110°, 120°, 130°, 150°
Material	SS316L, Titanium

NOTES

Blank lined area for notes.

C4 3.5mm Locking Plates

1. 3.5 mm Locking Compression Plates Low Contact
2. 3.5 mm Locking Reconstruction Plates Straight
3. 3.5 mm Locking Reconstruction Plates Curved
4. 3.5 mm Locking Reconstruction Plates Straight Low Profile
5. 3.5 mm Locking Reconstruction Plates Curved Low Profile
6. 3.5 mm Locking One Third Tubular Plates
7. 3.5 mm Locking T Plates 3 Head Holes
8. 3.5 mm Locking T Plates Angle 3 Head Holes
9. 3.5 mm Locking T Plates 4 Head Holes
10. 3.5 mm Locking Superior Anterior Clavicle Plates
11. 2.7/3.5 mm Locking Superior Anterior Clavicle Plates Lateral Extension
12. 3.5 mm Locking Clavicle Hook Plates
13. 3.5 mm Locking Medial Distal Tibia Plates with Tab
14. 3.5 mm Locking Medial Distal Tibia Plates without Tab
15. 3.5 mm Locking Proximal Tibia Plates
16. 3.5 mm Locking Medial Proximal Tibia Plates
17. 3.5 mm Locking Posteromedial Proximal Tibia Plates
18. 3.5 mm PHILOS Plates
19. 2.7/3.5 mm Locking Dorsolateral Distal Humerus Plates
20. 2.7/3.5 mm Locking Dorsolateral Distal Humerus Plates With Support
21. 2.7/3.5 mm Locking Medial Distal Humerus Plates
22. 3.5 mm Locking Extra Articular Distal Humerus Plates
23. 3.5 mm Locking Calcaneal Plates
24. 3.5 mm Locking Hook Plates
25. 3.5 mm Locking Olecranon Plates
26. 2.7/3.5 mm Locking Lateral Distal Fibula Plates
27. 3.5 mm Locking Anterolateral Distal Tibia Plates
28. 2.7/3.5 mm Locking Wrist Fusion Plates
29. 3.5 mm Locking Pediatric Hip Plates
30. 3.5 mm Locking Pilon Plates
31. 3.5 mm Locking Eight Plates

4.1 3.5 mm Locking Compression Plate Low Contact



Dia (mm)	3.5
Holes	6, 7, 8, 10, 12, 14, 16, 18, 20
Material	SS316L, Titanium

4.2 3.5 mm Locking Reconstruction Plate Straight



Dia (mm)	2.7 3.5 5
Holes	4, 6, 8, 10, 12, 14, 16, 18, 20
Material	SS316L, Titanium

4.3 3.5 mm Locking Reconstruction Plate Curved



Dia (mm)	2.7 3.5 5
Holes	6, 8, 10, 12, 14, 16, 18, 20
Material	SS316L, Titanium

4.4 3.5 mm Locking Reconstruction Plate Straight Low Profile



Dia (mm)	2.7 3.5 5
Holes	4, 6, 8, 10, 12, 14, 16, 18, 20
Material	SS316L, Titanium

4.5 3.5 mm Locking Reconstruction Plate Curved Low Profile



Dia (mm)	2.7 3.5 5
Holes	6, 8, 10, 12, 14, 16, 18, 20
Material	SS316L, Titanium

4.6 3.5 mm Locking One Third Tubular Plate



Dia (mm)	2.7 3.5
Holes	4, 6, 8, 10, 12, 14, 16, 18, 20
Material	SS316L, Titanium

4.7 3.5 mm Locking T Plate 3 Head Holes



Dia (mm)	3.5
Holes	3, 4, 5, 6
Material	SS316L, Titanium

4.8 3.5 mm Locking T Plate Angle 3 Head Holes



Size (mm)	3.5
Holes	3, 4, 5, 6
Direction	Left, Right
Material	SS316L, Titanium

4.9 3.5 mm Locking T Plate 4 Head Holes



Dia (mm)	3.5
Holes	3, 4, 5, 6, 7, 8
Material	SS316L, Titanium

4.10 3.5 mm Locking Superior Anterior Clavicle Plate



Size (mm)	3.5
Holes	6, 7, 8
Direction	Left, Right
Material	SS316L, Titanium

4.11 2.7/3.5 mm Locking Superior Anterior Clavicle Plate Lateral Extension



Size (mm)	3.5
Holes	4, 5, 6, 7, 8
Direction	Left, Right
Material	SS316L, Titanium

4.12 3.5 mm Locking Clavicle Hook Plate

Holes	4, 5, 6, 7
Hook Depth (mm)	12, 15, 18
Direction	Left, Right
Material	SS316L, Titanium

4.13 3.5 mm Locking Medial Distal Tibia Plate with Tab

Size (mm)	3.5
Holes	4, 6, 8, 10, 12, 14
Direction	Left, Right
Material	SS316L, Titanium

4.14 3.5 mm Locking Medial Distal Tibia Plate without Tab

Size (mm)	3.5
Holes	4, 6, 8, 10, 12, 14
Direction	Left, Right
Material	SS316L, Titanium

4.15 3.5 mm Locking Proximal Tibia Plate



Size (mm)	3.5
Holes	4, 6, 8, 10, 12, 14
Direction	Left, Right
Material	SS316L, Titanium

4.16 3.5 mm Locking Medial Proximal Tibia Plate



Size (mm)	3.5
Holes	4, 6, 8, 10, 12, 14, 16, 18
Direction	Left, Right
Material	SS316L, Titanium

4.17 3.5 mm Locking Posteromedial Proximal Tibia Plate



Size (mm)	3.5
Holes	2, 4, 6, 8, 10
Material	SS316L, Titanium

4.18 3.5 mm PHILOS Plate

Size (mm)	3.5
Holes	3, 4, 5, 6, 7, 8, 9, 10
Material	SS316L, Titanium

4.19 2.7/3.5 mm Locking Dorsolateral Distal Humerus Plates

Size (mm)	3.5
Holes	3, 5, 7, 9, 11, 13
Direction	Left, Right
Material	SS316L, Titanium

4.20 2.7/3.5 mm Locking Dorsolateral Distal Humerus Plates With Support

Size (mm)	3.5
Holes	3, 5, 7, 9, 11, 13
Direction	Left, Right
Material	SS316L, Titanium

4.21 2.7/3.5 mm Locking Medial Distal Humerus Plates



Size (mm)	3.5
Holes	3, 5, 7, 9, 11, 13
Direction	Left, Right
Material	SS316L, Titanium

4.22 3.5 mm Locking Extra-Articular Distal Humerus Plate



Size (mm)	3.5
Holes	4, 6, 8, 10, 12, 14
Direction	Left, Right
Material	SS316L, Titanium

4.23 3.5 mm Locking Calcaneal Plate



Size (mm)	3.5
Length (mm)	69, 76
Direction	Left, Right
Material	SS316L, Titanium

4.24 3.5 mm Locking Hook Plate



Size (mm)	3.5
Holes	3
Material	SS316L, Titanium

4.25 3.5 mm Locking Olecranon Plate



Size (mm)	3.5
Holes	2, 4, 6, 8, 10, 12
Direction	Left, Right
Material	SS316L, Titanium

4.26 2.7/3.5 mm Locking Lateral Distal Fibula Plate



Size (mm)	3.5
Holes	3, 4, 5, 6, 7
Direction	Left, Right
Material	SS316L, Titanium

4.27 3.5 mm Locking Anterolateral Distal Tibia Plate



Size (mm)	3.5
Holes	5, 7, 9, 11, 13, 15, 17
Direction	Left, Right
Material	SS316L, Titanium

4.28 2.7/3.5 mm Locking Wrist Fusion Plate




Size (mm)	3.5
Holes	8
Bend	Standard bend, Short bend, Straight
Material	SS316L, Titanium

4.29 3.5 mm Locking Pediatric Hip Plate



Size (mm)	3.5
Holes	100° (3 holes) 110° (3 holes) 120° (4 holes) 130° (3, 5, 7 holes) 150° (3 holes)
Angle	100°, 110°, 120°, 130°, 150°
Material	SS316L, Titanium

4.30 3.5 mm Locking Pilon Plate



Size (mm)	3.5
Holes	7, 9
Material	SS316L, Titanium

4.31 Eight Plate



Size (mm)	3.5
Length (mm)	12, 16, 20
Material	Titanium

NOTES

Blank lined area for notes.

C5 4.5/5mm Locking Plates

1. 4.5/5 mm Locking Narrow Compression Plate Low Contact
2. 4.5/5 mm Locking Broad Compression Plate Low Contact
3. 4.5/5 mm Locking Reconstruction Plate Straight
4. 4.5/5 mm Locking Reconstruction Plate Curved
5. 4.5/5 mm Locking T Plate 2 Head Holes
6. 4.5/5 mm Locking T Buttress Plate
7. 4.5/5 mm Locking Proximal Tibia Plate
8. 4.5/5 mm Locking Proximal Lateral Tibia Plate
9. 4.5/5 mm Locking Medial Proximal Tibia Plate
10. 4.5/5 mm Locking Distal Femoral Plate
11. 4.5/5 mm Locking Proximal Femoral Plate
12. 4.5/5 mm Locking Proximal Femoral Hook Plate
13. 4.5/5 mm Locking Condylar Femoral Plates
14. 4.5/5 mm Locking DHS Plate
15. 4.5/5 mm Locking DCS Plate 95 Degree
16. 4.5/5 mm Locking Trochanter Stabilizing Plate DHS
17. 4.5/5 mm Locking Paediatric Hip Plate

5.1 4.5/5 mm Locking Narrow Compression Plate Low Contact



Dia (mm)	4.5 / 5
Holes	6, 7, 8, 10, 12, 14, 16, 18, 20
Material	SS316L, Titanium

5.2 4.5/5 mm Locking Broad Compression Plate Low Contact



Dia (mm)	4.5 / 5
Holes	6, 7, 8, 10, 12, 14, 16, 18, 20
Material	SS316L, Titanium

5.3 4.5/5 mm Locking Reconstruction Plate Straight



Dia (mm)	4.5 / 5
Holes	4, 6, 8, 10, 12, 14, 16, 18, 20
Material	SS316L, Titanium

5.4 4.5/5 mm Locking Reconstruction Plate Curved



Dia (mm)	4.5 / 5
Holes	4, 6, 8, 10, 12, 14, 16, 18, 20
Material	SS316L, Titanium

5.5 4.5/5 mm Locking T Plate



Dia (mm)	4.5 / 5
Holes	3, 4, 5, 6, 8, 10
Material	SS316L, Titanium

5.6 4.5/5 mm Locking T Buttress Plate



Dia (mm)	4.5 / 5
Holes	3, 4, 5, 6, 8, 10
Material	SS316L, Titanium

5.7 4.5/5 mm Locking Proximal Tibia Plate



Dia (mm)	4.5 / 5
Holes	4, 6, 8, 10, 12, 14
Direction	Left, Right
Material	SS316L, Titanium

5.8 4.5/5 mm Locking Proximal Lateral Tibia Plate




Dia (mm)	4.5 / 5
Holes	5, 7, 9, 11, 13
Direction	Left, Right
Material	SS316L, Titanium

5.9 4.5/5 mm Locking Medial Proximal Tibia Plate



Dia (mm)	4.5 / 5
Holes	4, 6, 8, 10, 12, 14
Direction	Left, Right
Material	SS316L, Titanium

5.10 4.5/5 mm Locking Distal Femoral Plate



Dia (mm)	4.5 / 5
Holes	5, 7, 9, 11, 13
Direction	Left, Right
Material	SS316L, Titanium

5.11 4.5/5 mm Locking Proximal Femoral Plate



Dia (mm)	4.5 / 5
Holes	2, 4, 6, 8, 10, 12, 14, 16
Material	SS316L, Titanium

5.12 4.5/5 mm Locking Proximal Femoral Hook Plate



Dia (mm)	4.5 / 5
Holes	2, 4, 6, 8, 10, 12, 14, 16, 18
Material	SS316L, Titanium

5.13 4.5/5 mm Locking Condylar Femoral Plates



Dia (mm)	4.5 / 5
Holes	2, 4, 6, 8, 10, 12, 14, 16, 18
Direction	Left, Right
Material	SS316L, Titanium

5.14 5 mm Locking DHS Plate



DHS Plate Long Barrel (38 mm) Sizes

Angle	130°/ 135°/ 140°/ 145°
Holes	3 to 20
Material	SS316L, Titanium

DHS Plate Short Barrel (25 mm) Sizes

Angle	130°/ 135°/ 140°/ 145°
Holes	3 to 20
Material	SS316L, Titanium

5.15 Locking DCS Plate 95 Degre



Angle	95°
Holes	6 to 20
Material	SS316L, Titanium

5.16 Locking Trochanter Stabilizing Plate for DHS



Length (mm)	130
Material	SS316L, Titanium

5.17 5 mm Locking Paediatric Hip Plate



Size (mm)	5
Holes	100° (3 holes) 110° (3 holes) 120° (3 holes) 135° (3, 8, 7 holes)
Angle	100°, 110°, 120°, 135°
Material	SS316L, Titanium

NOTES

Blank lined area for notes.

C6 DHS-DCS & Angle Plates

1. Locking DHS Plate
2. DHS Plate
3. Locking Trochanter Stabilizing Plate for DHS
4. Locking DCS Plate 95 Degree
5. DCS Plate 95 Degree
6. Condylar Plate 95 Degree
7. Angled Blade Plates 130 Degree

6.1 Locking DHS Plate



DHS Plate Long Barrel (38 mm) Sizes

Angle	130°/ 135°/ 140°/ 145°
Holes	3 to 20
Material	SS316L, Titanium

DHS Plate Short Barrel (25 mm) Sizes

Angle	130°/ 135°/ 140°/ 145°
Holes	3 to 20
Material	SS316L, Titanium

6.2 DHS Plate



DHS Plate Long Barrel (38 mm) Sizes

Angle	130°/ 135°/ 140°/ 145°
Holes	3 to 20
Material	SS316L, Titanium

DHS Plate Short Barrel (25 mm) Sizes

Angle	130°/ 135°/ 140°/ 145°
Holes	3 to 20
Material	SS316L, Titanium

6.3 Locking Trochanter Stabilizing Plate for DHS



Length (mm)	130
Material	SS316L, Titanium

6.4 Locking DCS Plate 95 Degree



Angle	95°
Holes	6 to 20
Material	SS316L, Titanium

6.5 DCS Plate 95 Degree



Angle	95°
Holes	6 to 20
Material	SS316L, Titanium

6.6 Condylar Plate 95 Degree



Angle	95°
Blade Length (mm)	50/ 60/ 70/ 80
Holes	5 to 20
Material	SS316L, Titanium

6.7 Angled Blade Plate 130 Degree



Angle	130°
Blade Length (mm)	60/ 70/ 80/ 90/ 100
Holes	4 to 12
Material	SS316L, Titanium

C7 Bone Screws

1. 2 / 2.4 / 2.7 / 3.5 / 5 mm Locking Cortical Screws
2. 2 / 2.4 / 2.7 / 3.5 / 4.5 mm Cortical Screws
3. 3.5 / 4 / 5 / 6.5 mm Locking Cancellous Screws
4. 3.5 / 4 / 5 / 6.5 mm Cancellous Screws
5. 4 / 5 / 6.5 / 7.3 mm Locking Cannulated Screws
6. 4 / 4.5 / 6.5 / 7 / 7.3 mm Cannulated Screws
7. 2.5 / 3 / 3.5 / 4 / 4.5 / 5.5 / 6.5 / 7.5 Headless Compression Screws, Partially Thread
8. 2.5 / 3 / 3.5 / 4 / 4.5 / 5.5 / 6.5 / 7.5 Headless Compression Screws, Full Thread
9. Interlocking Nail Screws
10. Interference Screws
11. Herbert Screws Cannulated
12. Malleolar Screws
13. DHS-DCS Screws
14. Cross Head Screws
15. Anchor Screws

2 / 2.4 / 2.7 / 3.5 / 5 mm Locking Cortical Screws

available in **SS316L** & **Titanium**

any length can be manufacture

2 mm Locking Cortical Screws



2.4 mm Locking Cortical Screws



2.7 mm Locking Cortical Screws



3.5 mm Locking Cortical Screws



5 mm Locking Cortical Screws



2 / 2.4 / 2.7 / 3.5 / 4.5 mm Cortical Screws

available in **SS316L** & **Titanium**

any length can be manufacture

2 mm Locking Cortical Screws



2.4 mm Locking Cortical Screws



2.7 mm Locking Cortical Screws



3.5 mm Locking Cortical Screws



5 mm Locking Cortical Screws



3.5 / 4 / 5 / 6.5 mm Locking Cancellous Screws

available in **SS316L** & **Titanium**

any length can be manufacture

3.5 mm Locking Cancellous Screws



4 mm Locking Cancellous Screws



5 mm Locking Cancellous Screws, Part Thread



5 mm Locking Cancellous Screws, Full Thread



6.5 mm Locking Cancellous Screws, 16 mm Thread



6.5 mm Locking Cancellous Screws, 32 mm Thread



6.5 mm Locking Cancellous Screws, Full Thread



3.5 / 4 / 5 / 6.5 mm Cancellous Screws

available in **SS316L** & **Titanium**

any length can be manufacture

3.5 mm Cancellous Screws



4 mm Cancellous Screws, Full Thread



4 mm Cancellous Screws, Short Thread



6.5 mm Cancellous Screws, 16 mm Thread



6.5 mm Cancellous Screws, 32 mm Thread



6.5 mm Cancellous Screws, Full Thread



4 / 5 / 6.5 / 7.3 mm Locking Cannulated Screws

available in **SS316L** & **Titanium**

any length can be manufacture

4 mm Locking Cancellous Screws



4 mm Locking Cancellous Screws



6.5 mm Locking Cancellous Screws, 16 mm Thread



6.5 mm Locking Cancellous Screws, 32 mm Thread



6.5 mm Locking Cancellous Screws, Full Thread



7.3 mm Locking Cancellous Screws, 16 mm Thread



7.3 mm Locking Cancellous Screws, 32 mm Thread




7.3 mm Locking Cancellous Screws, Full Thread



4 / 4.5 / 6.5 / 7 / 7.3 mm Cannulated Screws

available in **SS316L & Titanium**

any length can be manufacture




4 mm Cancellous Screws, Short Thread

4 mm Cancellous Screws, Long Thread



4.5 mm Cancellous Screws, Short Thread


4.5 mm Cancellous Screws, Long Thread



6.5 mm Cancellous Screws, 16 mm Thread

6.5 mm Cancellous Screws, 32 mm Thread


6.5 mm Cancellous Screws, Full Thread



7 mm Cancellous
Screws, 16 mm Thread

7 mm Cancellous
Screws, 32 mm Thread

7 mm Cancellous
Screws, Full Thread



7.3 mm Cancellous
Screws, 16 mm Thread

7.3 mm Cancellous
Screws, 32 mm Thread

7.3 mm Cancellous
Screws, Full Thread

2.5 / 3 / 3.5 / 4 / 4.5 / 5.5 / 6.5 / 7.5 Headless Compression Screws, Partially Thread

available in **SS316L** & **Titanium**

any length can be manufacture

2.5 mm Headless Compression Screws, Partially Thread



3 mm Headless Compression Screws, Partially Thread



3.5 mm Headless Compression Screws, Partially Thread



4 mm Headless Compression Screws, Partially Thread



4.5 mm Headless Compression Screws, Partially Thread



5.5 mm Headless Compression Screws, Partially Thread



6.5 mm Headless Compression Screws, Partially Thread



7.5 mm Headless Compression Screws, Partially Thread



2.5 / 3 / 3.5 / 4 / 4.5 / 5.5 / 6.5 / 7.5 Headless Compression Screws, Full Thread

available in **SS316L** & **Titanium**

any length can be manufacture

2.5 mm Headless Compression Screws, Partially Thread



3 mm Headless Compression Screws, Partially Threaded



3.5 mm Headless Compression Screws, Partially Threaded



4 mm Headless Compression Screws, Partially Threaded



4.5 mm Headless Compression Screws, Partially Threaded



5.5 mm Headless Compression Screws, Partially Threaded



6.5 mm Headless Compression Screws, Partially Threaded



7.5 mm Headless Compression Screws, Partially Threaded



PFNA Blades



Length (mm)

70 / 75 / 80 / 85 / 90 / 95 / 100 / 105 / 110 / 115 / 120
125 / 130 / 135 / 140 / 145 / 150

available in **SS316L** & **Titanium**

PFNA2 Blades



Length (mm)

70 / 75 / 80 / 85 / 90 / 95 / 100 / 105 / 110 / 115 / 120
125 / 130 / 135 / 140 / 145 / 150

available in **SS316L** & **Titanium**

8 mm Proximal Cannulated Bolt



Length (mm)

70 / 75 / 80 / 85 / 90 / 95 / 100 / 105 / 110 / 115 / 120
125 / 130 / 135 / 140 / 145 / 150

available in **SS316L** & **Titanium**

6.4 mm Proximal Cannulated Bolt



Length (mm)

50 / 55 / 60 / 65 / 70 / 75 / 80 / 85 / 90 / 95 / 100 / 105 / 110
115 / 120 / 125 / 130 / 135 / 140 / 145 / 150

available in **SS316L** & **Titanium**

4.9 mm Locking Bolts



Length (mm)

22 / 24 / 26 / 28 / 30 / 32 / 34 / 36 / 38 / 40 / 42 / 44 / 46 / 48 / 50 / 52
54 / 56 / 58 / 60 / 65 / 75 / 80 / 85 / 90 / 95 / 100 / 105 / 110 / 115 / 120

available in **SS316L** & **Titanium**

3.9 mm Locking Bolt



Length (mm)

22 / 24 / 26 / 28 / 30 / 32 / 34 / 36 / 38 / 40 / 42 / 44 / 46 / 48 / 50 / 52
54 / 56 / 58 / 60 / 65 / 75 / 80 / 85 / 90 / 95 / 100

available in **SS316L** & **Titanium**

3.4 mm Locking Bolt



Length (mm)

22 / 24 / 26 / 28 / 30 / 32 / 34 / 36 / 38 / 40 / 42 / 44 / 46 / 48 / 50 / 52
54 / 56 / 58 / 60 / 65 / 75 / 80 / 85 / 90 / 95 / 100

available in **SS316L** & **Titanium**

Interference Screws Headless



Dia (mm)	6, 7, 8, 9, 10, 11
Length (mm)	20, 25, 30, 38
Material	Titanium, PEEK

Interference Screws Head

Dia (mm)	6, 7, 8, 9, 10, 11
Length (mm)	20, 25, 30, 38
Material	Titanium, PEEK



2.5 / 3 / 3.5 / 4 / 4.5 / 5.5 / 6.5 / Herbert Screws Cannulated

available in **SS316L & Titanium**

any length can be manufacture



3 mm Herbert Screws



3.5 mm Herbert Screws



4.5 mm Herbert Screws



5.5 mm Herbert Screws



6.5 mm Herbert Screws



3.5 / 4.5 mm Malleolar Screws

available in **SS316L** & **Titanium**
any length can be manufacture

3.5 mm Malleolar Screws



4.5 mm Malleolar Screws

**DHS-DCS Screws**

available in **SS316L** & **Titanium**
any length can be manufacture



Length (mm)

50, 55, 60, 65, 70, 75, 80, 85, 90, 95, 100, 105, 110, 115, 120, 125, 130, 135, 140, 145 and 150 mm

Cross Head Screws - Locking

2.0 mm Locking	6 mm - 30 mm
2.4 mm Locking	6 mm - 30 mm
2.7 mm Locking	6 mm - 40 mm

Cross Head Screws - Non Locking

1.5 mm Non Locking	6 mm - 30 mm
2.0 mm Non Locking	6 mm - 30 mm
2.4 mm Non Locking	6 mm - 30 mm
2.7 mm Non Locking	6 mm - 40 mm

Anchor Screws

available in **Titanium**

Dia (mm) 3 and 5

NOTES

Blank lined area for notes.

C8 CranioMaxilloFacial

1. Y Plate
2. Double Y Plate
3. X Plate
4. H Plate
5. Frame Plate Rectangle
6. Frame Plate Square
7. T Plate
8. Strut Plate
9. Strut Plate Curved
10. L Plate
11. L Plate Angle
12. Mesh Plate
13. TMJ Condylar Plate
14. TMJ Condylar Reconstruction Plate
15. Orbital Plate
16. Angle Plate Mandible
17. Crescent Plate Mandible
18. Straight Plate center space Mandible
19. Adaption Plate Mandible
20. Angle Reconstruction Plate Mandible
21. Double Angle Reconstruction Plate Mandible

8.01 Y Plate



Size (mm)	1.5 mm Non Locking 2 mm Locking 2.4 mm Locking 2.7 mm Locking
Material	Titanium

8.02 Double Y Plate



Size (mm)	1.5 mm Non Locking 2 mm Locking 2.4 mm Locking 2.7 mm Locking
Material	Titanium

8.03 X Plate



Size (mm)	1.5 mm Non Locking 2 mm Locking 2.4 mm Locking 2.7 mm Locking
Material	Titanium

8.04 H Plate



Size (mm)	1.5 mm Non Locking 2 mm Locking 2.4 mm Locking 2.7 mm Locking
Material	Titanium

8.05 Frame Plate Rectangle




Size (mm)	1.5 mm Non Locking 2 mm Locking 2.4 mm Locking 2.7 mm Locking
Material	Titanium

8.06 Frame Plate Square




Size (mm)	1.5 mm Non Locking 2 mm Locking 2.4 mm Locking 2.7 mm Locking
Material	Titanium

8.07 T Plate



Size (mm)	1.5 mm Non Locking 2 mm Locking 2.4 mm Locking 2.7 mm Locking
Holes	3 + 4 3 + 6
Material	Titanium

8.08 Strut Plate




Size (mm)	1.5 mm Non Locking 2 mm Locking 2.4 mm Locking 2.7 mm Locking
Holes	2 x 3 2 x 4
Material	Titanium

8.09 Strut Plate Curved




Size (mm)	1.5 mm Non Locking 2 mm Locking 2.4 mm Locking 2.7 mm Locking
Material	Titanium

8.10 L Plate



Size (mm)	1.5 mm Non Locking 2 mm Locking 2.4 mm Locking 2.7 mm Locking
Holes	2 +3 3+4 4+6
Direction	Left, Right
Material	Titanium

8.11 L Plate Angle



Size (mm)	1.5 mm Non Locking 2 mm Locking 2.4 mm Locking 2.7 mm Locking
Holes	2 +3 3+4 4+6
Direction	Left, Right
Material	Titanium

8.12 Mesh Plate



Thickness : 0.6 mm

Size :

50 x 50 mm

75 x 75 mm


100 x 100 mm

8.13 TMJ Condylar Plate




Size (mm)	1.5 mm Non Locking
	2 mm Locking
	2.4 mm Locking
	2.7 mm Locking
Material	Titanium

8.14 TMJ Condylar Reconstruction Plate



Size (mm)	1.5 mm Non Locking
	2 mm Locking
	2.4 mm Locking
	2.7 mm Locking
Material	Titanium

8.15 Orbital Plate




Size (mm)	1.5 mm Non Locking
	2 mm Locking
	2.4 mm Locking
	2.7 mm Locking
Holes	8 10 12
Material	Titanium

8.16 Angle Plate Mandible



Size (mm)	1.5 mm Non Locking
	2 mm Locking
	2.4 mm Locking
	2.7 mm Locking
Holes	3 + 3 4 + 4
Material	Titanium

8.17 Crescent Plate Mandible



Size (mm)	1.5 mm Non Locking
	2 mm Locking
	2.4 mm Locking
	2.7 mm Locking
Holes	2 + 2 3 + 3 5 7
Material	Titanium


8.18 Straight centre space Plate Mandible



Size (mm)	1.5 mm Non Locking
	2 mm Locking
	2.4 mm Locking
	2.7 mm Locking
Holes	2 + 2 3 + 3
Material	Titanium


8.19 Adaption Plate Mandible

Size (mm)	1.5 mm Non Locking 2 mm Locking 2.4 mm Locking 2.7 mm Locking
Holes	6 12 20 24 30
Material	Titanium




8.20 Angle Reconstruction Plate Mandible

Size (mm)	1.5 mm Non Locking 2 mm Locking 2.4 mm Locking 2.7 mm Locking
Holes	12+4 16+4 20+4 24+4
Material	Titanium



8.21 Double Angle Reconstruction Plate Mandible

Size (mm)	1.5 mm Non Locking 2 mm Locking 2.4 mm Locking 2.7 mm Locking
Holes	12+4 16+4 20+4 24+4
Material	Titanium




NOTES

Blank lined area for notes.

C9 Interlocking Nails

1. PFNA2 Nail / PFNA-II Nail
2. PFNA Nail
3. PFN Nail (Proximal Femoral Nail)
4. Expert Tibia Nail
5. Tibial (Tibia) Nail
6. Femur (Femoral) Nail
7. Antegrade Femoral Nail (AFN)
8. Femoral Reconstruction Nail (FRN)
9. Humeral Nail
10. TENS Nail (Titanium Elastic Nail)
11. Ender Nail
12. Rush Nail

available in **SS316L** & **Titanium**



PFNA2, Short Nail Sizes

Dia (mm)	Angle	Length (mm)
9 / 10 / 11 / 12	130°	180 / 200 / 240

PFNA2, Long Nail Sizes, Left & Right

Dia (mm)	Angle	Length (mm)
9 / 10 / 11 / 12 / 13 / 14	130°	260 / 280 /300/ 320/ 340/ 360/ 380/ 400/ 420

available in **SS316L** & **Titanium**




PFNA, Short Nail Sizes

Dia (mm)	Angle	Length (mm)
9 / 10 / 11 / 12	130°	180 / 200 / 240

PFNA, Long Nail Sizes, Left & Right

Dia (mm)	Angle	Length (mm)
9 / 10 / 11 / 12 / 13 / 14	130°	300 / 320 / 340 / 360 / 380 / 400 / 420

available in **SS316L** & **Titanium**




Proximal Femoral Nail (PFN), Short Nail Sizes

Dia (mm)	Angle	Length (mm)
9 / 10 / 11 / 12 / 13	130° / 135°	180 / 200 / 220 / 240

Proximal Femoral Nail (PFN), Long Nail Sizes, Left & Right


Dia (mm)	Angle	Length (mm)
9 / 10 / 11 / 12 / 13	130° / 135°	280 / 300/ 320/ 340/ 360/ 380/ 400/ 420

available in **SS316L** & **Titanium**



Femoral Recon Nail (Left & Right) Sizes

Dia (mm)	Length (mm)
9 / 10 / 11 / 12 / 13	280 / 300 / 320 / 340 / 360 / 380 / 400 / 420 / 440 / 460



available in **SS316L** & **Titanium**




Expert Tibia Nail Left & Right Sizes

Dia (mm)	Length (mm)
8 / 9 / 10 / 11 / 12 / 13	240/ 260/ 280 /300/ 320/ 340/ 360/ 380/ 400/ 420/ 440/ 460

4.9 mm Locking Bolt

avaliabile in **SS316L** & **Titanium**




Tibia Nail Left & Right Sizes

Dia (mm)	Length (mm)
8 / 9 / 10 / 11 / 12 / 13	240 / 260 / 280 / 300/ 320/ 340/ 360/ 380/ 400/ 420




available in **SS316L** & **Titanium**



4.9 mm Locking Bolt

Femur Nail Sizes

Dia (mm)	Length (mm)
8 / 9 / 10 / 11 / 12 / 13	280 / 300/ 320/ 340/ 360/ 380/ 400/ 420/ 440/ 460




4.9 mm Locking Bolt

available in **SS316L** & **Titanium**




Antegrade Femoral Nails (Left & Right) Sizes

Dia (mm)	Length (mm)
9 / 10 / 11 / 12 / 13 / 14	280 / 300/ 320/ 340/ 360/ 380/ 400/ 420/ 440/ 460




available in **SS316L** & **Titanium**



Humeral Nails Sizes

Dia (mm)	Length (mm)
7 / 8 / 9	140 / 160 / 180 / 200 / 220 / 240 / 260 / 280 / 300 / 320 /




Titanium Elastic Nail Sizes

Dia (mm)	Length (mm)
1.5	440
2.0	440
2.5	440
3.0	440
3.5	440
4.0	440
4.5	440
5	440


Ender Nail Sizes available in **SS316L** & **Titanium**

Dia (mm)	Length (mm)
2	240 / 260 / 280 / 300 / 320 / 340 / 360 / 380 / 400 / 420
2.5	240 / 260 / 280 / 300 / 320 / 340 / 360 / 380 / 400 / 420
3	240 / 260 / 280 / 300 / 320 / 340 / 360 / 380 / 400 / 420
3.5	240 / 260 / 280 / 300 / 320 / 340 / 360 / 380 / 400 / 420
4	240 / 260 / 280 / 300 / 320 / 340 / 360 / 380 / 400 / 420
4.5	240 / 260 / 280 / 300 / 320 / 340 / 360 / 380 / 400 / 420
5	240 / 260 / 280 / 300 / 320 / 340 / 360 / 380 / 400 / 420



Rush Nails Sizes available in **SS316L** & **Titanium**

Dia (mm)	Length (mm)
2	100 mm to 420 mm (10 mm difference)
2.5	100 mm to 420 mm (10 mm difference)
3	100 mm to 420 mm (10 mm difference)
3.5	100 mm to 420 mm (10 mm difference)
4	100 mm to 420 mm (10 mm difference)
4.5	100 mm to 420 mm (10 mm difference)
5	100 mm to 420 mm (10 mm difference)



NOTES

Blank lined area for notes.

C10 Hip Prosthesis

1. Bipolar Hip Prosthesis, Fenestrated
2. Bipolar Hip Prosthesis, Non Fenestrated
3. Austin Moore Hip Prosthesis
4. Thompson Hip Prosthesis
5. Lazcano Hip Prosthesis

10.1 Bipolar Hip Prosthesis, Fenestrated



Head Dia (mm)	37 to 57
Material	SS 316L

10.2 Bipolar Hip Prosthesis, Non Fenestrated



Head Dia (mm)	37 to 57
Material	SS 316L

10.1 Austin Moore Hip Prosthesis

Head Dia (mm)	37 to 57
Material	SS 316L



10.2 Thompson Hip Prosthesis

Head Dia (mm)	37 to 57
Material	SS 316L



10.2 Lazcano Hip Prosthesis

Head Dia (mm)	37 to 57
Material	SS 316L



NOTES

Blank lined area for notes.

C11 Spine System

1. Anterior Cervical Plate (ACLP)
2. 3.5 mm Anterior Cervical Screw
3. 4 mm Anterior Cervical Screw
4. Monoaxial Pedicle Screws
5. Monoaxial Reduction Screws
6. Polyaxial Pedicle Screws
7. Polyaxial Reduction Screws
8. Crosslink Connector / Transverse Connector Type 1
9. Crosslink Connector / Transverse Connector Type 2
10. Laminar Hook
11. Pedicle Hook
12. Connecting Rod
13. Cervical Cage | Anterior Cervical Interbody Fusion (ACIF Cage)
14. Lumbar Cage | Posterior Lumbar Interbody Fusion (PLIF Cage)
15. TLIF Cage | Transforaminal Lumbar Interbody Fusion Cage
16. Titanium Mesh Cage

11.1 Anterior Cervical Plate



Dia (mm)	Length (mm)
4	23
4	25
4	27
4	29
4	31
4	33
4	35
6	37
6	39
6	41
6	43
6	45
6	47
6	49
6	51
6	53
6	55
8	57
8	60
8	63
8	66
8	69
8	72
8	75
8	80

11.2 3.5 mm Anterior Cervical Screw



Length (mm) 12 / 14 / 16 / 18 / 20
available in **Titanium**

11.3 4 mm Anterior Cervical Screw



Length (mm) 12 / 14 / 16 / 18 / 20
available in **Titanium**

11.4 Monoaxial Pedicle Screw



Dia (mm) 4.5 / 5 / 5.5 / 6.5 / 7

Length (mm) 25 / 30 / 35 / 40 / 45 / 50 / 55 / 60

available in **Titanium**

11.5 Monoaxial Reduction Pedicle Screw



Dia (mm) 4.5 / 5 / 5.5 / 6.5 / 7

Length (mm) 25 / 30 / 35 / 40 / 45 / 50 / 55 / 60

available in **Titanium**

11.6 Polyaxial Pedicle Screws



Dia (mm) 4.5 / 5 / 5.5 / 6.5 / 7

Length (mm) 25 / 30 / 35 / 40 / 45 / 50 / 55 / 60

available in **Titanium**

11.7 Polyaxial Reduction Pedicle Screw



Dia (mm) 4.5 / 5 / 5.5 / 6.5 / 7

Length (mm) 25 / 30 / 35 / 40 / 45 / 50 / 55 / 60

available in **Titanium**

11.8 Crosslink Connector / Transverse Connector Type 1



Length (mm) 35 - 42 | 42 - 52 | 52 - 62

available in **Titanium**

11.9 Crosslink Connector / Transverse Connector Type 2



Length (mm) 80 / 90 / 100 / 120 / 150 / 200

available in **Titanium**

11.10 Connecting Rod



Dia (mm) 5.5 Length (mm) 80 to 400 (20 mm increment)

available in **Titanium**

11.11 Cervical Cage | Anterior Cervical Interbody Fusion (ACIF Cage)



Dia (mm) 4 / 5 / 6 / 7 / 8 / 9 / 10

Length (mm) 13

Width (mm) 16

available in **PEEK & Titanium**

11.12 Lumbar Cage | Posterior Lumbar Interbody Fusion (PLIF Cage)



Dia (mm) 8 / 10 / 12 / 14 / 16

Length (mm) 20 / 24 / 28 / 32

Width (mm) 10

available in **PEEK & Titanium**

11.13 TLIF Cage | Transforaminal Lumbar Interbody Fusion Cage



Dia (mm) 8 / 10 / 12 / 14 / 16

Length (mm) 24 / 28 / 32

available in **PEEK & Titanium**

11.14 Spine Mesh



Dia (mm) 10 / 12 / 14 / 16 / 18 / 20 / 22 / 24 / 26

Length (mm) 60 / 70 / 80 / 90 / 100

available in **Titanium**

C12 Pins / Staples / Wires

1. Steinmann Pins
2. Steinmann Pins Positive Profile Thread
3. Steinman Pin Centre Thread (Denham Pin)
4. Drill Bits
5. Schanz Screw (Schanz Pin)
6. Schanz Screw (Schanz Pin) Self Tapping
7. Schanz Screw (Schanz Pin) Tapered Thread
8. Kirschner Wire Single Trocar
9. Kirschner Wire Double Trocar
10. Kirschner Wire Double Trocar Fully Threaded
11. Suture Wire / Cerclage Wire
12. Cerclage Wire - Loop
13. Guide Wire for Tibia/Femur, Ball Tip, 2.5 x 900 mm
14. Guide Wire for Tibia/Femur, Plain, 2.5 x 900 mm
15. Gigli Saw Wire
16. Guide Wire for 4 mm Cancellous Cannulated Screws, 1.2 x 230 mm
17. Guide Wire for 6.5 mm Cancellous Cannulated Screws, 2 x 230 mm
18. Guide Wire for PFN 2/2.5/3 x 400 mm
19. Ilizarov Wire Bayonet Point
20. Ilizarov Wire Trocar Point
21. Ilizarov Olive Wire Bayonet Point
22. Ilizarov Olive Wire Trocar Point
23. Epiphyseal Staple
24. Coventry Staple

Steinmann Pins

Steinmann Pins - Single Tocar



Steinmann Pins - Single Thread



Steinmann Pins - Thread & Tocar



Dia (mm)	3 3.5 4 4.5 5 5.5 6 7 8
Length (mm)	100 to 250
Material	SS316L

Steinman Pin Positive Profile



Shaft Dia/ Thread Dia (mm)	4 - 5 4.5 - 5.5 5 - 6 5.5 - 6.5 6 - 7 6.5 - 8
Length (mm)	100 to 250
Material	SS316L

Steinman Pin Centre Thread (Denham Pin)



Dia (mm)	3 3.5 4 4.5 5 5.5 6 7 8
Length (mm)	100 to 250
Material	SS316L

Drill Bits

Cannulated Drill Bit



Cannulated Drill Bit, Quick Copuling



Drill Bit



Drill Bit, Quick Copuling



Dia (mm)	1 mm to 5 mm
Length (mm)	60 to 250
Material	SS316L

Schanz Screw / Schanz Pin Self Tapping



Dia (mm)	Thread Length (mm)	Total Length (mm)
4	20	80
4	30	100
4	40	125
4	40	150
4	40	175
4.5	20	80
4.5	30	100
4.5	40	125
4.5	40	150
4.5	40	175
5.0	30	100
5.0	40	125
5.0	60	150
5.0	60	175
5.0	80	200
5.0	80	250

Dia (mm)	Thread Length (mm)	Total Length (mm)
5.5	30	100
5.5	40	125
5.5	60	150
5.5	60	175
5.5	80	200
5.5	80	250
6	30	100
6	40	125
6	60	150
6	60	175
6	80	200
6	80	250
6.5	80	200
6.5	80	225
6.5	80	250
6.5	90	275
6.5	100	300

Schanz Screw / Schanz Pin - Self Drilling



Dia (mm)	Thread Length (mm)	Total Length (mm)
4	20	80
4	30	100
4	40	125
4	40	150
4	40	175
4.5	20	80
4.5	30	100
4.5	40	125
4.5	40	150
4.5	40	175
5.0	30	100
5.0	40	125
5.0	60	150
5.0	60	175
5.0	80	200
5.0	80	250
5.5	30	100
5.5	40	125
5.5	60	150
5.5	60	175
5.5	80	200
5.5	80	250
6	30	100
6	40	125
6	60	150
6	60	175
6	80	200
6	80	250
6.5	80	200
6.5	80	225
6.5	80	250
6.5	90	275
6.5	100	300

Schanz Screw / Schanz Pin - Tapered Thread



Dia (mm)	Thread Length (mm)	Total Length (mm)
4.5 - 3.5	20	80
4.5 - 3.5	20	100
4.5 - 3.5	30	100
4.5 - 3.5	30	120
4.5 - 3.5	30	140
4.5 - 3.5	40	120
4.5 - 3.5	40	140
4.5 - 3.5	40	160
6 - 5	30	100
6 - 5	30	120
6 - 5	30	140
6 - 5	30	160
6 - 5	40	140
6 - 5	40	160
6 - 5	40	180
6 - 5	40	200
6 - 5	50	160
6 - 5	50	180
6 - 5	50	200
6 - 5	50	220
6 - 5	60	180
6 - 5	60	200
6 - 5	60	220
6 - 5	60	240
6 - 5	70	200
6 - 5	70	220
6 - 5	70	240
6 - 5	80	220
6 - 5	80	250

K wire
Single Tocar



K wire
Double Tocar



K wire
Full Thread



K Wires

Dia (mm)	0.8 1 1.2 1.4 1.5 1.6 1.8 2 2.2 2.3 2.5 3
Length (mm)	80 100 125 150 175 200 225 250 300
Material	SS316L

K Wires Full Thread

Dia (mm)	1.5 2 2.5
Length (mm)	150 200
Material	SS316L

Suture Wire / Cerclage Wire

Dia (mm)	16 SWG (1.60 mm)
	18 SWG (1.25 mm)
	20 SWG (0.90 mm)
	22 SWG (0.70 mm)
	24 SWG (0.55 mm)
	26 SWG (0.45 mm)
	28 SWG (0.35 mm)
	30 SWG (0.30 mm)
Material	SS316L



Cerclage Wire - Loop



Guide Wire for Tibia/Femur, Ball Tip




Dia (mm)	2.5
Length (mm)	900
Material	SS316L

Guide Wire for Tibia/Femur



Dia (mm)	2.5
Length (mm)	900
Material	SS316L

Gigli Saw Wire



Guide Wire for 4 mm Cancellous Cannulated Screws



Dia (mm)	1.2
Length (mm)	230
Material	SS316L

Guide Wire for 6.5 mm Cancellous Cannulated Screws



Dia (mm)	2
Length (mm)	230
Material	SS316L

Guide Wire for PFN Nail (Proximal Femur Nail)



Dia (mm)	2 2.5 3
Length (mm)	400
Material	SS316L

Ilizarov Wire Bayonet Point



Dia (mm)	1.5 1.8
Length (mm)	400
Material	SS316L

Ilizarov Wire Trocar Point



Dia (mm)	1.5 1.8
Length (mm)	400
Material	SS316L

Olive Wire Bayonet Point




Dia (mm)	1.5 1.8
Length (mm)	400
Material	SS316L

Olive Wire Trocar Point




Dia (mm)	1.5 1.8
Length (mm)	400
Material	SS316L

Epiphyseal Staple



Dia A	Length B	Length C
1.5 mm	16 mm	19 mm
1.5 mm	22 mm	19 mm
1.5 mm	26 mm	19 mm
2 mm	16 mm	19 mm
2 mm	22 mm	19 mm
2 mm	26 mm	19 mm
2.5 mm	16 mm	19 mm
2.5 mm	22 mm	19 mm
2.5 mm	26 mm	19 mm

Coventry Staple



Dia D	Length A	Length B	Length C
2.5 mm	25 mm	33 mm	5 mm
2.5 mm	25 mm	33 mm	10 mm
2.5 mm	25 mm	33 mm	15 mm

C13 External Fixators

TUBULAR EXTERNAL FIXATOR

Schanz Screw / Schanz Pin Self Tapping
Schanz Screw / Schanz Pin - Self Drilling
Schanz Screw / Schanz Pin - Tapered Thread
Single Pin Clamp - Superior
Double Pin Clamp - Superior
Open Single Pin Clamp
Tube to Tube Clamp
Connecting Clamp
Self Holding Clamp
Twin Adjustable Clamp Straight
Twin Adjustable Clamp Curved
Transverse Clamp
Compression Distraction Device
Connecting Rod Dia 5 mm, L 200 mm, SS
Connecting Rod Dia 5 mm, L 200 mm, CF
Tubular Rod, Dia 11mm, Stainless Steel
Tubular Rod, Dia 11mm, Carbon Fiber
Triple Trocar Set with Handle
Sleeve for Cancellous Schanz Pin
Schanz Pin Introducer
T Wrench/Spanner 11 mm

SMALL TUBULAR EXTERNAL FIXATOR

Small Single Pin Clamp
Small Connecting Clamp
Small Open Single Pin Clamp
Connecting Rod Dia 4 mm Carbon Fiber
Connecting Rod Dia 4 mm Stainless Steel
Double Sleeve Dia 2.5 mm / 4 mm
T Wrench/Spanner 8 mm

RAIL EXTERNAL FIXATOR

Rail
Central Clamp
End Clamp
T Clamp
Ball and Socket Clamp
T Ball and Socket Clamp
Swivel Clamp
Graduated Swivel Clamp
Dyna Ring Clamp
Template for End & Central Clamp
Template for Swivel Clamp
Screw Guide
Drill Guide
Sandwich Plate
Bolt for Sandwich Plate
CD Unit
Half Pin Introducer
Allen Key 6 mm

ILIZAROV EXTERNAL FIXATOR

Full Ring Stainless Steel
Full Ring Carbon Fiber
5/8 Ring Stainless Steel
5/8 Ring Carbon Fiber
Half Ring Stainless Steel
Half Ring Carbon Fiber
3/4 Ring Carbon Fiber
1/4 Ring Carbon Fiber
Italian Femoral Arch 90 Degree, SS
Italian Femoral Arch 90 Degree Carbon Fiber
Italian Femoral Arch 120 Degree, SS
Italian Femoral Arch 120 Degree Carbon Fiber
Short Connection Plate SS
Connecting Plate Carbon Fiber
Connection Plates with Threaded Ends
Curved Plate
Twisted Plate
Threaded Socket
Threaded Rods
Threaded Rods Slotted
Telescopic Rod
Graduated Telescopic Rod
Oblique Support
Post Male
Post Female
Hinge - Male
Hinge - Female
Hinge-90°
Hinge Joint
Wire Fixation Bolt Slotted
Wire Fixation Bolt Cannulated
Single Pin Fixation Bolt
Ilizarov Rotation device
Pin Blocks
Multi Pin Fixation Clamp
Wire Fixation Buckle
Dynamic Wire Tensioners
Corticotome

13.1 Schanz Screw / Schanz Pin Self Tapping



Dia (mm)	Thread Length (mm)	Total Length (mm)
4	20	80
4	30	100
4	40	125
4	40	150
4	40	175
4.5	20	80
4.5	30	100
4.5	40	125
4.5	40	150
4.5	40	175
5.0	30	100
5.0	40	125
5.0	60	150
5.0	60	175
5.0	80	200
5.0	80	250
5.5	30	100
5.5	40	125
5.5	60	150
5.5	60	175
5.5	80	200
5.5	80	250
6	30	100
6	40	125
6	60	150
6	60	175
6	80	200
6	80	250
6.5	80	200
6.5	80	225
6.5	80	250
6.5	90	275
6.5	100	300

13.2 Schanz Screw / Schanz Pin - Self Drilling




Dia (mm)	Thread Length (mm)	Total Length (mm)
4	20	80
4	30	100
4	40	125
4	40	150
4	40	175
4.5	20	80
4.5	30	100
4.5	40	125
4.5	40	150
4.5	40	175
5.0	30	100
5.0	40	125
5.0	60	150
5.0	60	175
5.0	80	200
5.0	80	250
5.5	30	100
5.5	40	125
5.5	60	150
5.5	60	175
5.5	80	200
5.5	80	250
6	30	100
6	40	125
6	60	150
6	60	175
6	80	200
6	80	250
6.5	80	200
6.5	80	225
6.5	80	250
6.5	90	275
6.5	100	300

13.3 Schanz Screw / Schanz Pin - Tapered Thread



Dia (mm)	Thread Length (mm)	Total Length (mm)
4.5 - 3.5	20	80
4.5 - 3.5	20	100
4.5 - 3.5	30	100
4.5 - 3.5	30	120
4.5 - 3.5	30	140
4.5 - 3.5	40	120
4.5 - 3.5	40	140
4.5 - 3.5	40	160
6 - 5	30	100
6 - 5	30	120
6 - 5	30	140
6 - 5	30	160
6 - 5	40	140
6 - 5	40	160
6 - 5	40	180
6 - 5	40	200
6 - 5	50	160
6 - 5	50	180
6 - 5	50	200
6 - 5	50	220
6 - 5	60	180
6 - 5	60	200
6 - 5	60	220
6 - 5	60	240
6 - 5	70	200
6 - 5	70	220
6 - 5	70	240
6 - 5	80	220
6 - 5	80	250

13.4 Single Pin Clamp - Superior



13.5 Double Pin Clamp - Superior



13.6 Open Single Pin Clamp



13.7 Tube to Tube Clamp



13.8 Connecting Clamp



13.9 Self Holding Clamp



13.10 Twin Adjustable Clamp Straight



13.11 Twin Adjustable Clamp Curved



13.12 Transverse Clamp



13.13 Compression Distraction Device



13.14 Connecting Rod Dia 5 mm, Length 200 mm, Stainless Steel



13.15 Connecting Rod Dia 5 mm, Length 200 mm, Carbon Fiber



13.16 Tubular Rod, Dia 11mm, Stainless Steel



Length (mm) : 100, 150, 200, 250, 300, 350, 400, 450, 500

13.17 Tubular Rod, Dia 11mm, Carbon Fiber



Length (mm) : 100, 150, 200, 250, 300, 350, 400, 450, 500

13.18 Triple Trocar Set with Handle



13.19 Sleeve for Cancellous Shanz Pin



13.20 Schanz Pin Introducer



13.21 T Wrench/Spanner 11 mm



13.22 Small Single Pin Clamp

Size (mm)
4 x 2.5 | 4 x 4



13.23 Small Connecting Clamp

Size (mm)
4 x 2.5 | 4 x 4



13.24 Small Open Single Pin Clamp

Size (mm)
4 x 2.5 | 4 x 4



13.25 Connecting Rod Dia 4 mm Carbon Fiber



Length (mm) : 100, 150, 200, 250, 300, 350

13.26 Connecting Rod Dia 4 mm Stainless Steel



Length (mm) : 100, 150, 200, 250, 300, 350

13.27 Double Sleeve Dia 2.5 mm / 4 mm



13.28 T Wrench/Spanner 8 mm



13.29 Rail



Length (mm) : 240, 300, 400

13.30 Central Clamp



13.31 End Clamp



13.32 T Clamp



13.33 Ball and Socket Clamp



13.34 T Ball and Socket Clamp



13.35 Swivel Clamp



13.36 Graduated Swivel Clamp



13.37 Dyna Ring Clamp



13.38 Template for End and Central Clamp



13.39 Template for Swivel Clamp



13.40 Screw Guide



Length (mm) : 32, 45, 60, 80, 100

13.41 Drill Guide



Size (mm) : 3.2 X 40 | 3.2 X 60 | 3.2 X 80 | 4.8 X 40 |
4.8 X 60 | 4.8 X 80 | 4.8 X 120

13.42 Sandwich Plate



13.43 Bolt for Sandwich Plate



13.44 CD Unit



Length (mm) : 40, 80

13.45 Half Pin Introducer



13.46 Allen Key 6 mm



13.47 Full Ring Stainless Steel



Length (mm) : 100, 120, 140, 160, 180, 200, 240

13.48 Full Ring Carbon Fiber



Length (mm) : 100, 120, 140, 160, 180, 200, 240

13.49 5/8 Ring Stainless Steel



Length (mm) : 100, 120, 140, 160, 180, 200, 240

13.50 5/8 Ring Carbon Fiber



Length (mm) : 100, 120, 140, 160, 180, 200, 240

13.51 Half Ring Stainless Steel



Length (mm) : 100, 120, 140, 160, 180, 200, 240

13.52 Half Ring Carbon Fiber



Length (mm) : 100, 120, 140, 160, 180, 200, 240

13.53 3/4 Ring Carbon Fiber



Length (mm) : 100, 120, 140, 160, 180, 200, 240

13.54 1/4 Ring Carbon Fiber



Length (mm) : 100, 120, 140, 160, 180, 200, 240

13.55 Italian Femoral Arch 90 Degree, SS



Length (mm) : 100, 120, 140, 160, 180, 200, 240

13.56 Italian Femoral Arch 90 Degree Carbon Fiber



Length (mm) : 100, 120, 140, 160, 180, 200, 240

13.57 Italian Femoral Arch 120 Degree, SS



Length (mm) : 100, 120, 140, 160, 180, 200, 240

13.58 Italian Femoral Arch 120 Degree Carbon Fiber



Length (mm) : 100, 120, 140, 160, 180, 200, 240

13.59 Short Connection Plate SS



Holes : 2 to 10

13.60 Connecting Plate Carbon Fiber



Holes : 2 to 10

13.61 Connection Plates with Threaded Ends



Holes : 5, 7, 9, 11

13.62 Curved Plate



Holes : 3, 4, 5, 6, 8

13.63 Twisted Plate



Holes : 2, 3, 4

13.64 Threaded Socket



Length (mm) : 20, 30, 40, 60

13.65 Threaded Rods



13.66 Threaded Rods Slotted



Length (mm) : 75, 100, 125, 150

13.67 Telescopic Rod



Length (mm) : 100, 150, 200, 250, 300

13.68 Graduated Telescopic Rod



Length (mm) : 100, 150

13.69 Oblique Support



13.70 Post Male



Holes : 2, 3, 4

13.71 Post Female



Holes : 2, 3, 4

13.72 Hinge - Male



13.73 Hinge - Female



13.74 Hinge-90°



13.75 Hinge Joint



13.76 Wire Fixation Bolt Slotted



13.77 Wire Fixation Bolt Cannulated



13.78 Single Pin Fixation Bolt



13.79 Ilizarov Rotation device



13.80 Pin Blocks




Holes : 1, 2, 3, 4

13.81 Multi Pin Fixation Clamp



13.82 Wire Fixation Buckle



13.83 Dynamic Wire Tensioners



13.84 Corticotome



Size (mm) : 5, 7, 10, 15